

Polymerase Chain Reaction (PCR)

Materials

- ⌘ Sterile ddH₂O
- ⌘ 10X PCR Buffer
- ⌘ MgCl₂, 50 mM
- ⌘ dNTPs, 2.5 mM
- ⌘ 5' and 3' Primers, 2.5 μM
- ⌘ DNA Taq Polymerase, 5 u/μl (Bioline, Cat # BIO-21042)

Methods

⌘ PCR Master Mix

1. Keep all reagents on ice as well as the tubes for PCR.
2. Make a PCR master mix by adding the following ingredients in the order shown. Add the Taq Polymerase just before using. Mix well by vortexing.

Solution	Volume (μl)	Stock concentration	Final concentration
ddH ₂ O	35.5		
dNTPs	2	2.5 mM	0.1 mM
Primer 1	2	2.5 μM	0.1 μM
Primer 2	2	2.5 μM	0.1 μM
MgCl ₂	1.5 (0.5-5)	50 mM	1.5 (0.5-5) mM
PCR Buffer	5	10X	1X
Taq Polymerase	0.5	5 u/μl	2.5 u/reaction
Total	50		

The MgCl₂ concentration in the final reaction mixture is usually between 0.5 to 5.0 mM, and the optimum concentration is determined empirically (typically between 1.0 - 1.5 mM) Mg²⁺ ions. Generally, low Mg²⁺ leads to low yields (or no yield) and high Mg²⁺ leads to accumulation of nonspecific products (mispriming).

3. Add 1-5 μl DNA to PCR tube then add 20-24 μl master mix so that the total volume is 25 μl. The amount of DNA needed will depend upon the concentration; 10-100 ng of genomic DNA or 1-10 ng of plasmid DNA is sufficient.
4. Briefly spin the tubes to make sure all liquid is at the bottom.

⌘ PCR Reaction

1. Set up a PCR template such as;

Step	Temperature (°C)	Duration (minutes)	Number of cycles
Initial Hold	94	3 minute	1
Denaturation	94	30-90 seconds	25-40
Annealing	50-60	30-90 seconds	
Extension	72	30-90 seconds	
2 nd Hold	72	10	1

- Primers are usually 18-30 bp long with similar T_m values.
- 30 cycles is typical, too many can result in a high amount of nonspecific background.
- The initial hold step may be as long as 5 minutes.
- The denaturation, annealing, and extension time will depend upon on the length of the PCR fragment.
- The annealing temperature is usually 5°C lower than the true T_m of the primers.

References

Dieffenbach CW, Dveksler GS (eds) PCR Primer: A Laboratory Manual. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York

Erlich HA (ed) (1989) PCR Technology: Principles and Applications for DNA Amplification. Stockton Press, New York

Innis MA, Gelfand DH, Sninsky JJ (eds) (1995) PCR Strategies. Academic Press, San Diego, California

Innis MA, Gelfand DH, Sninsky JJ, White TJ (eds) (1990) PCR Protocols: A Guide to Methods and Applications. Academic Press, San Diego, California

Newton CR, Graham A (1994) PCR. BIOS Scientific Publishers, Oxford

Watson JD, Gilman M, Witkowski J, Zoller M (1992) Recombinant DNA, second edition. W.H. Freeman and Company, New York