GAC[®] 2100 G/GI







Safety Notices	1
Introduction	3
Sample Testing	3
Printer Connection (Optional)	0
Computer Interface	4
Features	4
GAC® 2100 G/GI Specifications	5
Installation and Setup	7
Unit Positioning	7
Disengaging Shipping Brace	8
Connecting AC Power	9
Removing Drawer Bottom (Optional)	9
Using Language Modes	10
Changing Access Codes	10
Using the Communication Ports	10
Start-IIn & Familiarization	11
Startup Procedure	
Statup Flocedure	10
Adjusting Screen Contract	IZ
Defining Keybeard Key Eurotions	12
Monu Scroope	12
Alarm	15
	15
Measuring Moisture	17
Measuring Moisture Preparing to Test Moisture	 17 17
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method)	 17 17 18
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers	 17 17 18 20
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters	17 17 18 20 20
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options	17 17 18 20 20 21
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page	17 17 18 20 20 21 22
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page	17 17 18 20 20 21 22 23
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu	17 17 18 20 20 21 22 22 23
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures	17 17 18 20 20 21 22 22 23 23 24
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key)	17 17 18 20 20 21 22 22 23 24 24
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key)	17 17 18 20 20 21 22 22 23 23 24 24 26
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key)	17 17 18 20 20 21 22 23 23 24 24 26 28
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key)	17 17 18 20 20 21 22 23 23 24 24 24 24 26 28 28
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Data (4 Key) Obtaining Basic Calibration Data (4 Key)	17 17 18 20 21 22 23 23 24 24 24 24 24 28 28 28 31
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key)	17 17 18 20 20 21 22 22 23 23 24 24 24 26 28 28 28 28 31
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key)	17 17 18 20 20 21 22 22 23 23 23 24 24 26 28 28 28 28 31 31 32
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key)	17 17 18 20 20 21 22 23 23 23 24 24 26 28 28 28 31 32 33
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key)	17 17 18 20 20 21 22 23 23 24 24 24 24 24 24 28 28 31 32 33 36
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key) Changing Date, Time and Date Format (5 Key)	17 17 18 20 21 22 23 23 23 24 24 24 24 24 24 28 28 31 32 33 36 37
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key) Changing the Labels (6 Key)	17 17 18 20 21 22 23 23 24 24 24 24 24 24 28 28 31 31 32 33 36 37 39
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key) Changing Date, Time and Date Format (5 Key) Changing the Labels (6 Key) Setting Communications Port Parameters (7 Key)	17 17 18 20 21 22 23 23 24 24 24 24 24 24 24 28 28 31 31 32 33 36 37 39 40

Testing Unit Performance	49
Nonactive Keys	49
Viewing COM1 Port (2 Key)	50
Viewing Com2 Transmit Port (3 Key)	50
Viewing Display Characters (5 Key)	51
Viewing Diagnostic Mode (6 Key)	52
Using Test Network Mode (7 Key)	54
Viewing/Printing Parameters (8 key)	56
Operator Maintenance	59
Cell Cleaning - Daily	59
Fuse Replacement	59
Error Codes	61
Appendix A - Communication Ports	65
COM1 Port Options	65
COM1 Pin Assignment Details	
COM2 Port Options	
COM2 Pin Assignment Details	67

Safety Notices	1
Introduction	3
Sample Testing	3
Printer Connection (Optional)	0
Computer Interface	4
Features	4
GAC® 2100 G/GI Specifications	5
Installation and Setup	7
Unit Positioning	7
Disengaging Shipping Brace	8
Connecting AC Power	9
Removing Drawer Bottom (Optional)	9
Using Language Modes	10
Changing Access Codes	10
Using the Communication Ports	10
Start-IIn & Familiarization	11
Startup Procedure	
Statup Flocedure	10
Adjusting Screen Contract	IZ
Defining Keybeard Key Eurotions	12
Monu Scroope	12
Alarm	15
	15
Measuring Moisture	17
Measuring Moisture Preparing to Test Moisture	 17 17
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method)	 17 17 18
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers	 17 17 18 20
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters	17 17 18 20 20
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options	17 17 18 20 20 21
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page	17 17 18 20 20 21 22
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page	17 17 18 20 20 21 22 23
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu	17 17 18 20 20 21 22 22 23
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures	17 17 18 20 20 21 22 22 23 23 24
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key)	17 17 18 20 20 21 22 23 23 24 24
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key)	17 17 18 20 20 21 22 22 23 23 24 24 26
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key)	17 17 18 20 20 21 22 23 23 24 24 26 28
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key)	17 17 18 20 20 21 22 23 23 24 24 24 24 26 28 28
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Data (4 Key) Obtaining Basic Calibration Data (4 Key)	17 17 18 20 21 22 23 23 24 24 24 24 24 28 28 28 31
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key)	17 17 18 20 20 21 22 22 23 23 24 24 24 26 28 28 28 28 31
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key)	17 17 18 20 20 21 22 22 23 23 23 24 24 26 28 28 28 28 31 31 32
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key)	17 17 18 20 20 21 22 23 23 23 24 24 26 28 28 28 31 32 33
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key)	17 17 18 20 20 21 22 23 23 24 24 24 24 24 24 28 28 31 32 33 36
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key) Changing Date, Time and Date Format (5 Key)	17 17 18 20 21 22 23 23 23 24 24 24 24 24 24 28 28 31 32 33 36 37
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key) Changing the Labels (6 Key)	17 17 18 20 21 22 23 23 24 24 24 24 24 24 28 28 31 31 32 33 36 37 39
Measuring Moisture Preparing to Test Moisture Measuring Moisture (Basic Method) Using Grain Sample ID Numbers Entering Numbers & Alpha Characters Using ID Options Viewing the Next Page Grain Selection Using Select Grain Menu Calibrating Procedures Viewing Calibration Data (1 Key) Entering/Changing Calibration Values (2 Key) Printing Calibration Grain List (3 Key) Obtaining Basic Calibration Data (4 Key) Establishing Parameters Clearing the Cycle Counter (1 Key) Choosing ID Options (2 Key) Configuring Output Options (3 Key) Changing the Access Code (4 Key) Changing Date, Time and Date Format (5 Key) Changing the Labels (6 Key) Setting Communications Port Parameters (7 Key)	17 17 18 20 21 22 23 23 24 24 24 24 24 24 24 28 28 31 31 32 33 36 37 39 40

Testing Unit Performance	49
Nonactive Keys	49
Viewing COM1 Port (2 Key)	50
Viewing Com2 Transmit Port (3 Key)	50
Viewing Display Characters (5 Key)	51
Viewing Diagnostic Mode (6 Key)	52
Using Test Network Mode (7 Key)	54
Viewing/Printing Parameters (8 key)	56
Operator Maintenance	59
Cell Cleaning - Daily	59
Fuse Replacement	59
Error Codes	61
Appendix A - Communication Ports	65
COM1 Port Options	65
COM1 Pin Assignment Details	
COM2 Port Options	
COM2 Pin Assignment Details	67



SAFETY NOTICES

Safety notices are one of the primary ways to call attention to potential hazards.



This Safety Alert Symbol identifies important safety messages in this manual. When you see this symbol, carefully read the message that follows. Be alert to the possibility of personal injury or death.

WARNING

Use of the word WARNING indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

ACAUTION

Use of the word CAUTION with the Safety Alert Symbol indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

CAUTION

Use of the word CAUTION without the safety alert symbol indicates a potentially hazardous situation which, if not avoided, may result in equipment damage.





INTRODUCTION

The DICKEY-john Grain Analysis Computer GAC[®] 2100 G/GI quickly tests grain samples for moisture content and other profile information. The unit loads a sample, weighs and checks temperature, runs test, displays results, and unloads the sample. The test results calculate automatically to display the moisture content and, if enabled, the temperature and approximate test weight (density) of the grain sample. Also, the product name, grain calibration date and sample ID number display for identification. An additional screen of information includes readings for the conductance (D1), capacitance (D2), weight (D3) and temperature (D4) values along with the currently set date and time.

Figure 1

GAC[®] 2100 G/GI Unit



SAMPLE TESTING

Performing a grain sample analysis consists of six basic steps;

- 1. Select a grain for testing through keyboard entries. Up to 16 grain choices are available on two pages.
- 2. Pour grain sample into the hopper (top of the unit) and enter a sample ID (optional) number.
- 3. Depress the LOAD key. The sample automatically loads (drops) the grain into the test cell and a strike-off arm levels and removes the excess.





Print key

- 4. Wait briefly for the test results to process. The moisture reading and grain temperature display.
- 5. Print out a hard copy of the test results. With an optional printer connected, the date and time with sample identification number and all test information prints out. By depressing the PRINT key, additional copies can be printed.
- 6. Press the **Unload** key to empty the test cell contents. The test cell rotates 180° (inverts) to dump the grain into a sample drawer or into a larger container below the counter top. After the cell returns to the upright position, the unit is ready for the next test sample.

PRINTER CONNECTION (OPTIONAL)

An optional printer connects to the GAC[®] 2100 G/GI to print out test data results. The printout may include the facility name and address, current date and time, sample ID number, product name, grain calibration date, percent moisture, grain temperature, unit serial number and the D1 through D4 values. The data can automatically print at the end of each measurement cycle or be manually initiated by pressing the **Print** key.

COMPUTER INTERFACE

The GAC[®] 2100 G/GI can output to a personal computer (PC) using the COM1 or COM2 ports. When using the COM2 port, a null cable or null adaptor must be used. A null adaptor is available from DICKEY-john.

FEATURES

Features of the GAC[®] 2100 G/GI include:

- Memory for up to 16 different grain calibrations.
- Keyboard entry for updating grain calibrations as revisions become necessary.
- Alpha/numeric sample identification numbers with automatic sequential numbering for record keeping purposes.
- Automatic ranging by switching to predetermined calibrations when moisture is above or below normal limits.
- Two serial output ports available for printer and/or external computer.
- Self-check mode to assure proper continual operation.
- Error messages display when out-of-limits moisture, grain weight or grain temperature occur.
- Help screens to assist operator's questions.



Figure 2





Top View

GAC® 2100 G/GI SPECIFICATIONS

- Supply Voltage and Frequency Limits: 85 Vac to 264 Vac 48-440 Hz at 1 Amp Max.
- **Operating Temperature Range:** 5° C (41° F) to 45° C (113° F). If the unit temperature range is exceeded, an ERROR message displays.
- Grain Temperature Range: 0° C (32° F) to 50° C (122° F). If the grain temperature range is exceeded, an ERROR message displays.
- Unit Temperature and Grain Temperature Difference Limit: The maximum allowable temperature differential between empty cell and grain under test is 20° C (36° F). If exceeded, an ERROR message displays.
- Classes/Types of Grain or Seed: Refer to the latest Calibration Bulletin.
- Weight: 11.8 kg (26 lbs) Shipping Weight: 15.0 kg (33 lbs)





INSTALLATION AND SETUP

Before shipping, the GAC[®] 2100 G/GI is inspected and tested for mechanical and electrical defects. After unpacking, visually inspect for damage occurring during transit. Save all packing materials until inspection is complete. If damage is found, file a claim with the carrier immediately. and notify your DICKEY-john sales representative.

UNIT POSITIONING

Choose a clean environment for the GAC[®] 2100 G/GI that is protected from rapid changes in ambient temperature. Avoid a hazardous (classified) location as defined in Article 500 of the NFPA Handbook of the National Electrical Code.

Use the self-contained level on the panel top to position the GAC[®] 2100 G/GI on a flat level or nearly level surface (See Figure 2). Adjust the four rubber foot pads to level and stabilize the unit (See Figure 3). When moving the unit to a different location, re-leveling may be necessary. After adjustment is complete, all four leveling foot pads should be set to the lowest level possible with one pad fully seated against the bottom panel. This assures minimum height above the counter top.

Figure 3

Leveling Unit with Rubber Foot Pads



To adjust a foot pad height, use a flat wrench to unlock the upper nut from the bottom panel of the selected foot pad. Then, adjust the height as required.

After adjusting, use a Phillips screwdriver in the center of the footpad to hold the screw from turning and tighten the upper nut against the bottom panel of the unit. Ensure the lower nut remains tight against the footpad.



DISENGAGING SHIPPING BRACE

The shipping brace must be moved to the operate position before applying power. Since the grain test cell cycles through an unload sequence during power up, damage may occur with the shipping brace in the transport position.

The GAC[®] 2100 G/GI is shipped with the grain test cell supported by an internal shipping/transport brace. This brace is to support the test cell whenever the unit is transported. The position of the brace is controlled by a handle, located on the rear panel behind the grain test cell (refer to Figure 4). A small spool is placed over the handle during shipping to prevent possible bumping and releasing of the handle. Remove this spool and store as part of shipping materials. Push the handle in to unlock and rotate counterclockwise to the locked operate position as indicated on the back panel. The movement of the brace is viewable through the front window of the unit.



Viewing Unit from Rear



Rotate the brace back into supporting position before moving or shipping the unit. If the cell is not adequately supported, operating performance can be adversely effected due to cell bouncing. Some effort is required to bring the brace back into position. It may be helpful to reach inside the drawer front and lift up on the cell while attempting to lock the brace into the transport position. Place the spool over the handle to prevent the accidental release of the shipping brace.

CONNECTING AC POWER

Verify that the ON/OFF (I/O) switch is in the OFF (O) position and that the shipping brace is in the OPERATE position. Plug the line cord into the connector on the rear panel next to the ON/OFF switch. Plug the male end into an appropriate 3-wire (grounded) outlet.



The grounding pin on the line cord connects directly to the $GAC^{\textcircled{R}}$ 2100 G/GI frame. When using an adapter, make sure the grounding wire is connected properly to a good earth ground to prevent a shock hazard.

REMOVING DRAWER BOTTOM (OPTIONAL)

Normal use of the sample drawer can be bypassed by removing the drawer bottom and cutting a hole into the work counter to allow the tested grain samples to fall through to a larger container below. The size of the hole must be at least equal to the drawer dimension and located directly beneath the grain sample drawer.

Remove the bottom of the grain sample drawer by carefully lifting the two plastic tabs at the front of the drawer (refer to Figure 5). Using a small screwdriver, gently pry up on each tab and slide the bottom panel out the rear of the drawer. Be careful not to break the tabs.

Figure 5

Removing Drawer Bottom



When finished, slide the grain drawer back into position. This prevents grain from spilling onto the work surface in front of the unit.

If replacing the drawer bottom later, slide the bottom toward the drawer front until it latches under the two plastic tabs.



CHANGING ACCESS CODES

An access code is necessary to use certain screens. Whenever the code is required, a screen appears stating ENTER ACCESS CODE. The factory preset access code is set to "0". To change the access code for security reasons, refer to CHANGING THE ACCESS CODE in the SETTING UP PARAMETERS section.

USING THE COMMUNICATION PORTS

Two communication ports are located on the rear of the unit transfer data. These are: (a) Printer Output Port (COM1) and (b) Computer Output Port (COM2). Each must be properly configured to communicate with external devices. Port configuration is described in ESTABLISHING PARAMETERS.



START-UP & FAMILIARIZATION

This section defines the keyboard layout, display messages, and other operator controls.

START-UP PROCEDURE

The unit performs a series of self checks every time the unit is powered on to determine the status of the (1) load/strike-off mechanism, (2) unload function, (3) empty test cell weight, and (4) conductance/capacitance measurement circuitry. The self-checks begin when the hopper doors open, the strike-off arm swings across the cell and back again, and the cell dumps. When the cycle ends (about 15 seconds), the MAIN menu displays.

Before starting the procedure below, ensure the ON/OFF (I/O) switch is in the OFF (O) position and that the shipping brace is in the OPERATE position. These are both located on the rear of the unit (refer to Figure 4).

1. Turn the power switch ON. Initially, a language screen appears showing eight choices as shown in Figure 6.

Figure 6

Start-up Screens



2. Multiple languages are pre-programmed into the unit. There are 8 choices available on the start-up screen. Select a language by pressing the key number as indicated next to the language. This screen only remains available for approximately eight seconds. Once language is selected, an asterisk appears and remains in front of the selection until changed in another start-up cycle.



- 3. Wait for the remaining portion of the cycle to finish. A second screen will appear and display PLEASE WAIT while performing self checks. Software revision numbers will display on bottom two lines of screen (refer to Figure 6).
- 4. After the MAIN menu displays, proceed with the following paragraphs to preview general operating features.

DISPLAY SCREEN

The display screen is a graphics Liquid Crystal Display (LCD) and is backlit for low light level viewing.

ADJUSTING SCREEN CONTRAST

The Display Contrast Control is on the rear panel above the COM1 and COM2 ports (refer to Figure 4). Contrast Control changes the contrast between the displayed characters and the background. Adjust the control using a flat blade screwdriver as necessary for a comfortable contrast level in varying light conditions.

DEFINING KEYBOARD KEY FUNCTIONS

All Operator Controls are on the front panel except the Power ON/ OFF switch and Contrast Control which are located on the rear panel (refer to Figure 4).

The sealed membrane keyboard contains 20 pressure sensitive keys requiring light finger pressure for key actuation. Valid key closures result in short tone bursts while invalid keystrokes generate longer tone bursts. Each key function or key group is defined as follows and illustrated in Figure 7.

Numbers (1 through 0)

Decimal point (.)

Minus (-) Keys

All numeric keys perform dual functions: (1) make MENU selections and (2) enter numeric values. The four numeric keys with small arrows have the additional function of screen navigation. When making a numeric entry, a prompt on the display identifies the value entered and a cursor shows the destination of each entry. The minus (-) sign is accepted at any point in the entry and displays at the front of the entered string.



Figure 7

Keyboard Layout



Each entry immediately displays as the cursor moves to the right in preparation for the next entry. Only the first entry of the decimal point is acknowledged.



Backspace key



Next Page key

The **Backspace** key moves the cursor to the left to delete the last keystroke or to clear an entire entry field.

The **Next Page** key advances to the next page when the <more> prompt appears in the lower right corner of the display as shown in Figure 8. The NEXT PAGE key sequences the display forward to additional pages of information.

Figure 8

Additional Pages Screen



OPERATOR'S MANUAL





Previous Page key



Enter key



Load key



Help key



Print key



1^{*} [&] 3^{*}





Up and Down keys

The **Previous Page** key sequences backward through previously displayed pages. Neither the Previous Page or Next Page keys roll-over the pages to the beginning or end.

The Enter key accepts the character string after an entry is complete.

The **Load** key initiates a grain cell loading sequence to begin a test cycle.

The **Help** key provides context sensitive help screens. If the <more> message appears at the lower right of the screen, additional pages are accessible with the Next Page key. To use the HELP INDEX in the HELP screen, press the 0 key and then select the desired topic. To resume normal operation, press the Help key a second time.

The **Print** key initiates a printout of current information on the screen (i.e. grain moisture analysis or calibrations data). If the automatic printout option is selected, the Print key only initiates duplicate copies. Help information does not print. COM SETUP options allows the Print key to transmit to either the printer port or the communications/modem port. For example, the system may be configured to automatically print results through the printer port while the Print key transfers results to the computer port.

Pressing the **Unload** key instructs the test cell to dump its contents after a grain test is finished.

The **Left and Right** Keys move the cursor horizontally during alpha entries and certain functions during the SETUP mode.

The **Up and Down** Keys move the cursor vertically during alpha entries and certain functions in the SETUP mode.



MENU SCREENS

The MAIN menu appears immediately after all start-up checks successfully end. The MAIN menu is the home screen that all other screens are accessed. The MAIN menu displays five major selections or categories (refer to Figure 9).

A number will appear next to each MAIN menu selection and the message SELECT ITEM NUMBER will appear at display bottom. Select the desired number to access the appropriate menu. Additional menu levels will be available for some selections. Refer to Figure 10 for a general configuration of the Menu Tree layout.

Figure 9

MAIN menu Selections



ALARM

An audible alarm is enabled for 0.5 seconds at the end of each grain analysis or for an error condition. A shorter beep of 0.1 second sounds for each valid key stroke.



Figure 10

Menu Tree





MEASURING MOISTURE

The MAIN menu appears immediately after all start-up checks have successfully concluded (refer to START-UP & FAMILIARIZATION). The MAIN menu displays five major selections or categories. This section describes the MEASURE MOISTURE procedure.

PREPARING TO TEST MOISTURE

To measure moisture:

1. Choose MEASURE MOISTURE (1 Key) from the MAIN menu to display page one of the SELECT GRAIN menu as shown in Figure 11.

Figure 11







Next Page key



Previous Page key

2. Locate the grain type for testing and select the corresponding number to the left. Press the **Next Page** key if desired grain is not shown on first screen. After selecting a grain, the unit cycles through self-check tests to ensure that the program memory and calibrations are valid. Subsequent moisture measurements default to the previouslyMAIN menu selected grain. If a different grain is desired, press the **Previous Page** key to return to the MAIN menu. Press the 2 key, SELECT GRAIN, and make a new selection (refer to SELECTING GRAINS for additional information). NOTE: Pressing "NEXT PAGE" displays the



second set of eight grains that can be tested with the factory calibration setup.

3. Enter a sample ID number. This number identifies the sample tested on the paper printout. If the ENTER REQUIRED FOR ID is turned on (refer to CHOOSING ID OPTIONS in the ESTABLISHING PARAMETERS section), the ENTER SAMPLE ID screen displays and an ID number is required (refer to Figure 12). Other ID numbering options skip this screen. For now, enter any number and press the ENTER key (i.e. Press 1). Further explanation on using ID numbers are described later in this section.

Figure 12

Screen Sequence Using Enter Required for ID



MEASURING MOISTURE (BASIC METHOD)

After all selected grain calibration self-checks are finished (approximately 10 seconds), the message FILL SAMPLE HOPPER, ENTER ID NO. THEN...PRESS LOAD TO BEGIN will display (refer to Figure 13).

To measure moisture:

- 1. Pour grain into the sample hopper. Fill the hopper until the grain is at least 1.25 cm (1/2 inch) above the separation line between the hopper top and hopper compartment (heap grain slightly). The exact volume is unimportant except enough grain must be present to overfill the cell (approx .5 liter or 2 cups). Excess grain spills over the cell top and falls into the grain drawer below.
- 2. If an ID number has not already been entered, enter a number or skip to step 3.
- 3. Press the **Load** key on the keyboard to start the test cycle. After pressing the Load key, the hopper doors open to drop the grain into the test cell. The cell fills and the striker arm swings across the top of the test cell to wipe away excess grain.



Enter key



Load key



Do not insert fingers or other foreign objects into the opened hopper doors.

 Wait approximately 15 seconds for moisture test to complete as shown in Figure 13. A screen appears with message, MEASURING -PLEASE WAIT.

Figure 13



FILL SAMPLE HOPPER, ENTER ID NO. THE	
PRESS LOAD TO BEG	IN
WHEAT HRW 9 ID:AD1234	60918
Press LOAD Key	
	\rightarrow
MEASURING	
PLEASE WAIT	
WHEAT HRW 9	60918
ID: AD38	
13.4% ^{MOISTI}	URE
WHEAT HRW	960918
TEMP: 30°C DENSITY: 77.9 kg/hl ID: AD38	<more></more>
Press DUMP Key	
	\sim
PERFORMING SELF CHECKS	
PLEASE WAIT	
WHEAT HRW	960918



Next Page key



- 5. The grain name and percent moisture will appear on the display along with the sample temperature in degrees Centigrade (Fahrenheit), Density, Sample ID, and Grain Calibration issue date.
- 6. The option <more> is displayed in the lower right corner of the screen. Press the **Next Page** key to view additional pages. In this case, the D1 through D4 values display along with the present date (lower left corner) and time (lower right corner). Details of this screen are described later.
- 7. When measurement is complete, press the **Unload** key. The test cell inverts, dumping the grain and then returns to the upright position ready for the next measurement cycle.

 Start a new moisture measurement for the same grain type. The GAC[®] 2100 G/GI remembers the grain type until changed or power is turned off. The SELECT GRAIN menu does not appear again until selected.

USING GRAIN SAMPLE ID NUMBERS

Grain samples can be identified with three optional ID numbers:

- 1. Automatic number sequencing (i.e. 1, 2, 3...)
- 2. Manual entry (numbers can be manually entered or skipped)
- 3. ENTER REQUIRED FOR ID (produces an extra screen, ENTER SAMPLE ID) which prompts a required entry before proceeding (See Figure 12). The ID field can be either numeric or alphanumeric. If numeric only, the digits 0-9, the decimal point and the minus sign are entered with those keys on the keyboard. Other characters are entered from an alpha screen.

ENTERING NUMBERS & ALPHA CHARACTERS

Entering alpha characters (letters) into certain screens, such as an ENTER SAMPLE ID number, is accomplished through a special alpha screen. Press the NEXT PAGE key to access the alpha screen from the ENTER SAMPLE ID screen. The **Next Page** key only functions if an alpha entry is possible.

To view and/or use the alpha screen:

- Prepare for a moisture test. When the PRESS LOAD TO BEGIN message appears, observe the ID field in the lower left corner below the grain type. A small cursor is blinking requesting either a number or letter. Numbers are entered directly from the key pad and completed by pressing the Enter key.
- 2. To enter letters, press the **Next Page** key to access the alpha screen. The alpha screen shows the alphabet and special characters. In the lower left corner, the word TEXT is displayed. (refer to Figure 14).



Next Page key



Enter key



Backspace key

20 / MEASURING MOISTURE







Backspace key



Previous Page key

Figure 14

Alpha Screen

ABCDEFGHIJKLMNOPQRSTUVQXYZ _ #\$%&'() *,+
ÇŐŹÄĆÉĹĽŚÖÜŤŁAŇŽĘČÁÂĚS ĂĐĎËŇIÍTŮÓßOŃŚŔÚŰ
ݧŘ
TEXT: -

- 3. Move the vertical pointer (center of screen) to the first letter. The small vertical pointer identifies the selected character. The pointer is moved about the screen using the four arrow keys (left/right and up/down).
- 4. After selecting a letter, press the **Enter** key. The letter appears next to the word TEXT on the lower corner. Repeat the process until the phrase is complete. If an error is made, press the **Backspace** key and re-enter the corrected character.
- 5. After completing the text string, press **Previous Page** key. The text string now becomes part of the ID field on the PRESS LOAD TO BEGIN screen. Additional numbers may be inserted after the text string. Deleting is accomplished by entering the alpha screen again and using the **Backspace** key to delete letters.

USING ID OPTIONS

ID numbers are typically used to identify grain samples. The numbers can be a simple number (i.e. 1, 2, 3, etc.), a combination of letters and numbers (i.e. operator's initials plus sample number) or perhaps the customer's name (i.e. BROWN). The ID field may contain up to 16 alpha numeric characters.

The ID option is selected through the SETUP MENU, ID OPTIONS (refer to SETTING UP PARAMETERS). ID numbers are added to each tested sample using one of three methods:

- 1. **Inserted manually** Any ID field, including none at all, can be assigned to grain samples before each test. This flexibility allows users to develop individualized numbering schemes.
- Automatic Sequencing Numbers can be sequentially assigned automatically to each moisture test. For example, 1, 2, 3, etc. or specific letters (i.e. initials) can be placed in front of the numbers (i.e. JW12). Only the numbers at the end of the field change automatically.
- 3. Entry required for ID An ID field is required before each moisture test can proceed. This method ensures each sample is properly identified usually for printer output. A screen appears before a moisture test begins asking the user to ENTER SAMPLE ID. Alpha characters can be used if desired.





Next Page key

VIEWING THE NEXT PAGE

After the grain analysis finishes and the moisture content displays, the <more> appears in the lower right corner of the screen indicating another page is available. Press the **Next Page** key to view the calculated values (D1, D2, D3 and D4) determined from the grain measurement (refer to Figure 15). These values are relatively unimportant to the operator but are important to service personnel.

A general definition of each value is as follows;

- D1 Conductance measurement to indicate grain surface moisture.
- D2 Capacitance measurement to indicate internal grain moisture.
- D3 Weight measurement to indicate grain weight.
- D4 Temperature to indicate grain temperature.

Figure 15

Next Page Typical Calculated D-Values

$\langle \rangle$	\			
	D1: nnnn			γ
	D2: nnnn			
	D3: nnnn			
	D4: nnnn			
	TEST WEIGHT OFFSET:	+0.4 LB/BU		
	DATE: 05/18/96		TIME: 16:42:20	
	ID: 12			人
$\langle \rangle$				$\overline{)}$



GRAIN SELECTION

The second item on the MAIN menu is the SELECT GRAIN menu used to select specific grain types for testing. The SELECT GRAIN menu automatically appears for the first moisture test each time power is applied or can be selected anytime a different grain type is to be tested. After the first test, the SELECT GRAIN menu must be selected manually.

USING SELECT GRAIN MENU

The GAC[®] 2100 G/GI stores up to 16 different grain calibrations.

To select a grain for testing:

- 1. Choose the SELECT GRAIN menu from the MAIN menu by pressing the 2 key. Page 1 of the SELECT GRAIN menu shows the first eight grain calibrations available. The lower left corner of the the screen displays the page number and the lower right indicates more pages (refer to Figure 16).
- 2. Press the **Next Page** key to advance to the next page. The next eight grain calibration names appear. The <MORE> prompt disappears.
- 3. Press the PREVIOUS PAGE key to move backward. The **Previous Page** key is the counterpart to the **Next Page** key.
- 4. Press the number key corresponding to the number adjacent to the desired grain. After selecting the grain, the GAC[®] 2100 G/GI automatically enters the MEASURE MOISTURE mode and retains this grain selection until a new grain is selected again or the power is turned off.

Figure 16

Select Grain Menus Page 1



5. Press the **Next Page** key to view the second page of eight grains as shown in Figure 17.



Next Page key



Previous Page key



Figure 17

Select Grain Menus Page 2



CALIBRATING PROCEDURES

The CALIBRATION MENU allows the operator or authorized personnel to:

- · view the calibration constants of selected grains
- · enter new or change existing calibrations
- print the entire list of stored grain calibrations
- obtain calibration data necessary to generate new calibrations.

Selecting a function is achieved by pressing the number key corresponding to the desired function as shown in Figure 18.

Figure 18

Calibration Menu with Status Page



VIEWING CALIBRATION DATA (1 KEY)

The VIEW CALIBRATION menu allows viewing and printing of a selected grain calibration. (refer to Figure 19). Editing cannot be performed in this function.



Next Page key



Print key

To View Grain Calibrations:

- 1. From the CALIBRATION MENU, choose VIEW CALIBRATION by pressing the 1 key. Page 1 of 2 pages of grains will be displayed.
- 2. Choose a grain calibration to view from the grain list and select the corresponding number key. Use the **Next Page** key to view page 2 of the grain listings.
- 3. The VIEW CALIBRATION screen shows the Grain Name, Moisture Range, Calibration Issue Date and the nine Calibration K values.
- 4. If a printer is connected, press the **Print** key to print out a copy of this screen. If AUTO RANGE is specified in the calibration, this information is also printed. The lower right corner of the VIEW CALIBRATION screen indicates <more>meaning another page is available.
- 5. Press the Next Page key to display the AUTO RANGE screen. The AUTO RANGE screen shows the grain name and calibrations automatically used to extend the moisture range if the normal maximum and minimum limits are exceeded. Separate grain calibrations or grains can be referenced for each limit.

Figure 19

View Calibration Menu Showing Wheat HRW





ENTERING/CHANGING CALIBRATION VALUES (2 KEY)

To Enter/Change Calibration Values:

- 1. From the Calibration menu, press the 2 key to access the Enter/ Change Calibration scrren.
- Select the grain to be altered or locate any empty position to enter a new calibration. Select the appropriate key corresponding to the listing. A screen will appear displaying TO PROCEED, UNSEAL ACCESS PORT ON REAR PANEL AND PRESS WHITE PUSHBUTTON. Refer to Figure 20 for additional information.

IMPORTANT: Do not perform this action unless you are authorized

Figure 20

Changing Calibrations







Next Page key



Enter key



Previous Page key



Print key

- 3. If selecting an existing calibration, a large flashing block will appear on the grain selection. The screen will display the current K values for that grain. The large flashing block signifies that alpha characters may be entered from the alpha screen. Pressing the **Backspace** key, as prompted on the bottom of the screen, causes the bottom prompt to disappear and the large flashing block to change to a cursor (flashing line). When the cursor appears, the alpha screen may be accessed by pressing the **Next Page** key. If numeric entries are desired, those values are made directly with the number keys.
- 4. Use the Alpha screen to enter a grain name (up to 16 characters). Press the **Enter** key after each character selection to place the character on the text line.
- 5. Press the **Previous Page** key to return to the Enter/Change Calibration screen. The new grain title is now in place.
- 6. Press the **Enter** key to accept the new grain title. The cursor will advance to the next position (Issue Date).
- 7. Enter the Calibration Issue Date (up to 8 characters) and press the Enter key. The cursor will return to a flashing block in the ISSUE position. The block signifies the date may include alpha characters (i.e. 15May06), but normally is entered as numbers in a format such as YYMMDD (i.e. 060515). Pressing the Backspace key will cause the flashing block to disappear, allowing the Next Page key to be used for selecting the Alpha screen.
- 8. Press the **Enter** key to accept the date. The cursor will appear as a line (flashing) under the first digit of K1 (constant). The flashing cursor (line) signifies this entry accepts only numbers.
- Enter the K1 constant. When entering the first digit the remaining spaces are blank. All four digits of the constant must be entered when making a change. If more than four digits of the constant must be entered, a long beep will indicate an illegal entry.
- 10. When finished, press the **Enter** key to accept the constant and to advance to the next position. If the screen is accidentally advances before finishing, press the **Enter** key as many times as necessary to cycle through all menu items.
- 11. Enter K2 through K9 constants in the same manner as K1. Press the **Enter** key after entering each constant (4 digits).
- 12. Press the **Print** key to print the calibration. The present screen values will print.
- 13. Press the **Next Page** key to enter the Auto Range calibration selection. The Auto Range screen allows entry for grains using higher or lower moisture levels than the current calibration. These values are automatically selected if the grain moisture is outside the standard setting.
- 14. Press the **Backspace** key and choose a grain calibration to be used if the sample moisture is above the upper limit. After selecting a grain, the display will return to the Auto Range screen. Repeat this procedure again for the lower limit grain calibration.
- 15. Press the 0 (zero) key to delete either auto range calibration. Press the **Previous Page** key to return to the Grain Calibration screen.

PRINTING CALIBRATION GRAIN LIST (3 KEY)

A printout of all grain calibrations stored in memory can be obtained by pressing the **Print** key (See Figure 21).

To Print a Calibration List:

- 1. Press the 3 key from the CALIBRATION MENU. The PRINT CALIBRATION LIST screen will appear.
- 2. Press the **Print** key. The list of all grain calibrations will print.

Figure 21

Printing Calibration Lists





Next Page key



Enter key



Load key

OBTAINING BASIC CALIBRATION DATA (4 KEY)

This mode allows the user to generate data for developing grain calibrations. The samples used for calibration should span the desired moisture range.

In order to obtain good calibration data, the samples must be as evenly distributed throughout the moisture range as possible. The grain samples should extend at least 2% past each end of the desired moisture range, should contain no spoilage or condensation, and be free of broken pieces or foreign matter. The GAC[®] 2100 G/GI readings obtained may be sent to DICKEY-john to derive the calibration constants for the new grain.

To Obtain Accurate Calibration Constants:

- 1. Turn the GAC[®] 2100 G/GI power switch On and allow approximately thirty minutes to warm up before taking sample readings. Turn the printer On and confirm that the GAC[®] 2100 G/GI is setup and functioning properly.
- 2. From the MAIN menu, select the CALIBRATION MENU by pressing the 3 key. Then select OBTAIN CALIBRATION DATA by pressing the 4 key.
- 3. Enter a Sample Identification Number so the sample can be properly identified. Key in a sample ID number (use **Next Page** key if alpha ID desired). Press the **Enter** key.
- 4. Pour the sample into the hopper and press the Load key.



NOTE: The Facility Name and Location, Date/Time, GAC[®] 2100 G/GI Serial Number and CALIBRATION DATA only print when the first sample is run in this mode. Refer to Figure 23).

- 5. After the measurement cycle is complete, the screen displays four D values (refer to Figure 22). If using a printer, a printout automatically results. If a printer is not being used, record the D-values and sample ID.
- 6. Write the grain type under the CALIBRATION DATA. Write the actual moisture percentage of each sample (obtained by an accurate approved means) next to STD MOIST. Run additional samples.
- 7. Run as many samples as possible, making sure to write in the actual moisture percentage for each sample. After running each sample, press the UNLOAD key.

Figure 22





Figure 23

Calibration Data Printout

DICKE) 10653	/-JOHN
04/21/98	06:42:33
S/N: 1640)-10653
CALIBF	RATION DATA
ID:	1
STD. MOIS	ST:
	D1: 0 D2: 14 D3: 2 D4: 179



ESTABLISHING PARAMETERS

The SETUP MENU establishes basic operating parameters for the unit. These include resetting the cycle counter, ID options, output options, setting the access code, setting date and time, entering labels, and configuring communication ports.

From the MAIN menu, the SETUP MENU is accessed by pressing the 4 key (refer to Figure 24). From the SETUP MENU, eight sub-menus will appear. Most functions require the access code before changes can be made. If the access code has not been changed, the factory installed temporary code is set at 0 (zero).

Figure 24

Setup Menu Choices





Enter key

CLEARING THE CYCLE COUNTER (1 KEY)

The Cycle Counter increments once each time a moisture measurement or data collection sequence is successful. This screen displays the total counts accumulated since the last reset.

To reset the Cycle Counter:

- 1. From the SETUP MENU, press the 1 key as shown in Figure 25. A screen will appear displaying ENTER ACCESS CODE.
- 2. Enter the current personal access code (or 0) and press the **Enter** key. Enter 0, if the code has not been changed from the factory setting.
- 3. Press Enter again to reset the number (count) to zero. The counter will return to zero.

NOTE: The date, time, and cycle counter values are automatically sent to the COM1 and COM2 ports each time the counter is cleared, resulting in a printout of the most recently cleared counter contents. The printout should be retained for records.


Clearing the Cycle Counter



CHOOSING ID OPTIONS (2 KEY)

ID OPTIONS determine how each moisture measurement sample is identified on the output data. ID numbers appear in the lower left corner of each test result-screen. All sample ID numbers may contain up to 16 alphanumeric characters. Three methods of sample identification are available for display:

- AUTOMATIC SEQUENTIAL ID Automatically increments once for each measurement cycle but resets each time power is turned off.
- ENTER REQUIRED FOR ID With each measurement, a screen appears prompting ENTER SAMPLE ID before the grain sample can be measured. This option ensures the operator assigns a grain sample identification number/label for each test printout.
- OPTIONAL ENTRY ID If both the above selections are checked NO, the ID number can be added to each measurement manually or not at all. Manually added ID numbers do not require any specific format or sequencing.

To configure ID Options:

- From the SETUP MENU, press the 2 key to select the ID OPTIONS screen (refer to Figure 26). A screen will appear displaying ENTER ACCESS CODE.
- Enter your personal access code (or 0) and press Enter. The ID OPTION screen displays only two lines but allows three choices – The AUTOMATIC SEQUENTIAL ID and ENTER REQUIRED FOR ID functions cannot both be enabled (YES) at the same time.



Enter key





Up and Down keys



Select the desired function. Use the Up/Down Arrow keys (0 and 5 keys) to select and the Left/Right Arrow keys (1 and 3 key) to set the desired value. If both selections are set to NO, grain sample ID is optional by the operator for each test.

Figure 26

Selecting ID Options



CONFIGURING OUTPUT OPTIONS (3 KEY)

OUTPUT OPTIONS configure information appearing on the display/printer for the measured test results (refer to Figure 27). These selections are typically made once and not changed again but can be altered as necessary. These selections are only available to authorized personnel as outlined below.

To Change the Output Options:

- 1. From the SETUP MENU, select the OUTPUT OPTIONS (3 key). A screen will appear displaying instructions on how to proceed.
- 2. If access is authorized, locate and press the small white pushbutton under the sealed access panel at the rear of the unit. The panel is located in the lower left rear corner (viewed from rear) of the unit. The pushbutton is on the circuit board assembly inside the cavity.
- Verify each option and set as necessary. Three pages must be verified as shown in Figure 27. Use Up/Down Arrow keys (0 or 5) to select functions and Left/Right Arrow keys (1 or 3) to set each function to the appropriate setting.



Setting Output Options



Options include:

- **DISPLAY & OUTPUT TEST WEIGHT:** Enables/disables (YES/NO) approximate grain test weight to display and print.
- DISPLAY & OUTPUT TEMPERATURE: Enables/disables (YES/NO) grain temperature to display and print.

DICKEY-john®

- DISPLAY OUT-OF-LIMITS RESULT: Error messages flash on the display each time tested grain parameters exceed the normal range of that particular grain. If enabled, three standard parameters are displayed – moisture, grain temperature and approximate grain test weight. Each measured parameter value appears next to its respective legend on the display. The display out-of-limit options selected below establish whether or not the enabled values appear on the display. Eight separate settings are available for suppressing the enabled grain parameters.
- 1. **YES** Enabled readings display for all out-of-limit conditions.
- NO Enabled readings do not display for any out-of-limit condition.
 M Enabled readings display if moisture is out-of- limits but do not
- display if either weight and/or temperature is out-of-limits.
- 4. **W** Enabled readings display if weight is out-of-limits but do not display if either moisture and/or temperature is out-of-limits.
- 5. **T** Enabled readings display if temperature is out-of-limits but do not display if either weight and/or moisture is out-of-limits.
- 6. **MW** Enabled readings display if moisture and/or weight is out-of-limits but do not display if temperature is out of limits.
- 7. **MT** Enabled readings display if moisture and/or temperature is out-of-limits but do not display if weight is out-of-limits.
- 8. **WT** Enabled readings display if weight and/or temperature is out-of-limits but do not display if moisture is out-of-limits.
- **OUTPUT OUT-OF-LIMITS RESULT**: If enabled, three standard parameters are printed moisture, grain temperature and approximate grain test weight. Each measured parameter value appears next to its respective legend on the printout. The output out-of-limit options selected establish whether or not the enabled values print. These choices are identical to the display out-of-limits conditions above and are selected in a similar fashion.
- **OUTPUT REPORT OF ALL ERRORS:** Enables all internal checked errors to print.
- **MOISTURE PRINTOUT RESOLUTION**: Selects the resolution to either one tenth (0.1) percent or one hundredth (0.01) percent for moisture printout.
- **RADIX POINT:** Selects a decimal point (.) or comma (,) to display and printout.
- **TEMPERATURE FORMAT:** Selects either degrees C (°C) or degrees F (°F) to display and print out for temperature.
- WEIGHT FORMAT: Selects either pounds per bushel (lb/bu) or kilograms per hectoliter (kg/hl) to display and printout.



NOTE: The access code is set to 0 when shipped from the factory. After selecting a different code, record the new number, and store for safekeeping. DO NOT lose the new code.



Enter key

CHANGING THE ACCESS CODE (4 KEY)

The Access Code establishes a security level for nearly all Setup Menu functions.

To change the access code:

- 1. From the SETUP MENU, select the EDIT ACCESS CODE (4 key). The screen will appear displaying ENTER ACCESS CODE.
- 2. Enter the present code (or 0 if still set for factory code) and press the **Enter** key. The EDIT ACCESS CODE screen will display a blank (to be filled in) and the present Access Code displays at right (refer to Figure 28).
- 3. Enter the new Access Code. Choose any numeric code up to 4-digits. Enter a number and press **Enter**. The new number must be used from this point forward to enter the access code level secured screens.

Figure 28

Changing the Access Code

	SETUP 1: CLEAR COUNTER 2: ID OPTIONS 3: OUTPUT OPTIONS 4: EDIT ACCESS CODE SELECT ITEM NUMBER	MENU 5: DATE/TIME MENU 6: LABELS MENU 7: COM MENU 8: SERVICE MENU (NEXT PAGE: STATUS)
EM	EDIT ACCESS C	
	EDIT A	ICCESS CODE

CHANGING DATE, TIME AND DATE FORMAT (5 KEY)

This function establishes correct date and time for display and print out with each grain moisture measurement.

To Set the Time:

- 1. From the SETUP MENU, select DATE/TIME MENU (5 key). Three options are available as shown in Figure 29:
- SET TIME
- SET DATE
- SET DATE FORMAT

Figure 29

Setting Time



- 2. From the DATE/TIME MENU, select SET TIME using the 1 key. The SET TIME screen will appear and display ENTER ACCESS CODE.
- 3. Enter Access Code to gain access and press **Enter**. The time displays based upon a 24 hour clock.
- 4. Change the time by entering new numbers from the keyboard. The time format is in hours, minutes, seconds (HH.MM.SS). Each unit is two digits and separated by a decimal point. For example, to enter 3:32:59 P.M., enter 15.32.59. To separate use the decimal point. For accuracy, set the time slightly ahead and wait. At the exact time, press the **Enter** key to start the clock.





Enter key

To change the date:

- 1. From the DATE/TIME MENU, select SET DATE by pressing the 2 key as shown in Figure 30. The SET DATE screen will display ENTER ACCESS CODE.
- 2. Enter the access code to gain access to the SET DATE menu and press the Enter key. Enter numbers in the format shown at the bottom of the display; Month.Day.Year. Each date component is two digits in length and separated by a decimal point (i.e. to enter 5/15/06, enter 05.15.06). The number will appear on the bottom left of screen. Press the **Enter** key to accept the corrected date.

Figure 30







Left and Right keys



Previous page key

To Change the Date Format:

- 1. From the DATE/TIME MENU, select SET DATE FORMAT by pressing the 3 key (refer to Figure 31). The screen will display ENTER ACCESS CODE.
- 2. Enter Your Access Code to gain access to the SET DATE menu, then press the Enter key. Enter numbers in the format shown at the bottom of the display; Month.Day.Year. Each date component is two digits in length and separated by a decimal point (i.e. to enter 5/15/06, enter 05.15.06). The number appears on the bottom left of screen. Press the Enter key to accept the corrected date.
- Press the left or right arrow keys (1 or 3) to select one of three formats. When finished, press the **Previous Page** key to accept the new format and return to DATE/TIME MENU. Press the **Previous Page** key once more to return to the SETUP MENU.







CHANGING THE LABELS (6 KEY)

The LABELS MENU allows editing of User Name and attaches to the grain moisture measurements. The basic procedure is outlined in the following five steps. Specific details and illustrations for each label follow the basic outline (refer to Figure 31).

• USER NAME – appears on each measurement printout (not the display) and may be the organization's name. The User Name may be two lines and up to 16 characters each.



Enter key



Backspace key



Previous Page key

To Create Labels:

- 1. From the SETUP MENU, select the LABELS MENU (6 key).
- 2. Select the desired label by pressing the key number next to the label type. The ENTER ACCESS CODE screen will appear on display.
- 3. Enter the access code and press the **Enter** key. For the User Name screen (refer to Figure 32), a large flashing block highlights the first character to be edited. Press the **Backspace** key as prompted in the lower bottom corner of the display. The prompt disappears and is replaced with a cursor (underline) beneath the first character of the first line.
- 4. Enter the selected characters/numbers to reflect the label. Numbers and letters may be mixed by using the alpha screen.
- 5. When completed, press the **Previous Page** key and return to the LABELS MENU.



Editing User Name



SETTING COMMUNICATIONS PORT PARAMETERS (7 KEY)

Selection of the COM (Communications) MENU is accessed through the SETUP MENU by pressing the 7 key. From this menu, the operator selects output parameters for communicating with external devices such as a PC computer and/or a printer. Two port connectors appear on the rear of the unit directly above the sealed access cover. The COM1 (Data Communications Equipment) port connector is above the COM2 (Data Terminal Equipment) port connector. The COM MENU contains eight items with number 5 nonactive (refer to Figure 33). Selections 2, 3 and 4 configure the COM1 port while 6, 7 and 8 configure the COM2 port. Each of the COM MENU items is edited in a similar fashion.







To Configure the Communications Ports

(Basic entry procedure – Follow-up to specific subject for finishing details):

- 1. From the MAIN menu, press the 4 key (SETUP MENU) followed by the 7 key (COM MENU). The COMMUNICATIONS MENU will appear with seven active choices and one nonactive selection.
- 2. Choose the desired function by pressing the key number next to the choice. The ENTER ACCESS CODE screen will appear.
- 3. Enter the access code and press the **Enter** key. The selected COM screen will appear for editing.
- 4. Proceed with the necessary editing. The selections are:
 - COM PORT SELECT
 - COM1 or COM2 FORMAT
 - COM1 or COM2 SERIAL DATA SETUP
 - COM1 or COM2 CONTROL LEAD SETUP

To Configure COM PORT SELECT press the 1 key from the COMMUNICATIONS MENU.

1. After completing Steps 1 through 4 from the basic entry procedure, edit the configuration. Four items appear on the menu for editing and configuring output information of both ports (COM1 and COM2). Each item can be set to NONE, COM1, COM2 or COM1/COM2 (both ports) using the **Left/Right Arrow** keys to cycle through the choices. The pointer in front of each item can be moved with the **Up/Down Arrow** keys.







GAC2100 G/GI Grain Analysis Computer

11001-1381-200605







The four outputs are defined as follows (refer to Figure 34):

- AUTOMATIC RESULTS OUTPUT the configured port produces an output to automatically printout the moisture results.
- MANUAL RESULTS OUTPUT produces a printout for the configured port of the results by pressing the **Print** key. Repeatedly pressing the PRINT key results in duplicate copies.
- ERROR REPORT OUTPUT the configured port produces a printout of errors occurring during internal self checks. Out-of-limit errors of measured results do not print.
- SETUP AND TEST OUTPUT the configured port produces a printout of Setup and Test data by pressing the Print key.
- 2. Once editing is complete, press the **Previous Page** key to return to the COMMUNICATIONS MENU and for the next selection.

To Configure COM1 or COM2 FORMAT press 2 or 6 key from the COMMUNICATIONS MENU):

 After completing Steps 1 through 4 from the Basic Entry Procedure above, edit the configuration. Both the COM1 FORMAT and COM2 FORMAT selections are similar and are described together. The COM FORMAT SELECTION Menu displays two lines each containing several choices as shown in Figure 35.



Print key



Previous Page key





Up & Down Arrow keys



Left & Right Arrow keys



Previous Page key

- 2. Select the OUTPUT FORMAT using the pointer (refer to Figure 35). Place the pointer in front of the correct line with the **Up/Down Arrow** keys.
- 3. Edit the line by using the left/right arrow keys to cycle through six variations. Select the appropriate format from the following six choices.
 - PRN20 Standard 20 column printer output.
 - PRN80 Standard 80 column printer output.
 - ST20 20 column test output of printer parameters.
 - PRN81 Standard 80 column printer output (2 lines).
 - CSV01 Comma separated variables.
 - PRN21 Special 20 column printer output. PRN 20 format with the Grain Calibration Date Code, Instrument Serial Number and D-values removed.
- 4. Select the LINE TERMINATOR with the pointer. An RS-232 line of input data must be terminated with one of the following characters or character pairs.
 - LF (Line Feed)
 0A Hex
 - CR (Carriage Return) 0D Hex
 - LF + CR 0A, 0D Hex
 - CR + LF
 0D, 0A Hex





5. After editing is complete, press the **Previous Page** key to return to the COMMUNICATIONS MENU for the next selection.





Previous Page key

To Configure COM1 or COM2 SERIAL DATA SETUP (3 and 7 keys from the Communications Menu):

Both the COM1 and COM2 SERIAL DATA SETUP selections are similar and therefore described together.

- After completing Steps 1 through 4 from the Basic Entry Procedure above, edit the configuration. Either of the COM SERIAL DATA SETUP Menus display four lines but only the top line (Baud rate) can be changed for the COM1 SERIAL DATA SETUP (refer to Figure 36). Select the baud rate line for the COM1 SERIAL DATA SETUP of the printer. Use the Left/Right Arrow keys (1 or 3 key) to cycle the Baud Rate through seven choices; 300, 600, 1200, 2400, 4800, 9600, and 19200 bits/ second. Choose the one that matches the rate of the equipment in use.
- 2. Select the DATA BITS, PARITY and STOP BITS in that order. These items for the COM1 port cannot be changed. The COM2 port items can be adjusted. Definitions of each item is as follows:
 - DATA BITS Set for 7 or 8 bit character lengths.
 - PARITY Set to NONE, EVEN or ODD via left/right arrow keys.
 - STOP BITS Set to 1 or 2.

Figure 36

Accessing the COM1 Serial Data Setup Screen



3. After editing is complete, press the **Previous Page** key to return to the COMMUNICATIONS MENU for the next selection.

DICKEY-john



Previous Page key

To Configure COM1 or COM2 CONTROL LEAD SETUP (4 or 8 Keys from the CONTROL LEAD COMMUNICATIONS MENU):

Both the COM1 and COM2 CONTROL LEAD SETUP selections are similar and therefore described together in this procedure (See Figure 37). Differences are noted as applicable.

 After completing Steps 1 through 4 from the Basic Entry Procedure above, edit the configuration of the COM1 Port screen. The COM1 (DCE) CONTROL LEAD SETUP Menu displays two lines (See Figure 36, Accessing COM CONTROL LEAD SETUP Screens). The Busy

Line, BUSY (11), notifies the GAC[®] 2100 G/GI when the buffer in the printer (computer) is full and cannot accept more data. The polarity of this line is important and is set to be compatible with the printer. The Control Lead Busy Sense has three polarity choices: IGNORE, ACTIVE (+) or ACTIVE (-). The DTR (20) (Data Terminal Ready) line notifies the GAC[®] 2100 G/GI that the printer (computer) is turned on. The Control Lead DTR Sense has two states: IGNORE or ACTIVE.

- Edit the configuration of the COM2 Port screen. The COM2 (DTE) CONTROL LEAD SETUP Menu displays seven items (See Figure 36, Accessing COM CONTROL LEAD SETUP Screens) but only two lines can be changed.
 - RTS (4) : ON (+) or AUTO
 - DSR (6) : IGNORE or ACTIVE
- 3. After editing is complete, press the **Previous Page** key to return to the COMMUNICATIONS MENU for the next selection.







Figure 38

Common Acronym Definitions

- BUSY Line busy
- CTS Clear To Send
- DCE Data Communications Equipment
- DSR Data Set Ready
- DTE Data Terminal Equipment

- DTR Data Terminal Ready
 - Ring Indicator

RI

- RLSD Receive Line Signal Detect
- RTS Request To Send
- RxD Receive Data
- TxD Transmit Data

DICKEY-john®

Nonactive (5 Key from COMMUNICATIONS MENU)

The nonactive selection (refer to Figure 33) has no function and only emits a long beep when pressed. No other action results.

USING THE SERVICE MENU (8 KEY)

The SERVICE MENU contains individual unit calibration constants and is accessed by DICKEY-john trained and equipped service personnel.





TESTING UNIT PERFORMANCE

The TESTS MENU displays diagnostic tests and status information for operating values. The information is useful to technicians for overall performance evaluations and isolating problems. During normal operation, these are not generally referenced. However, problem identification can be assisted by reviewing certain values occasionally and observing changes. Most values cannot be edited but the contents can be printed for records. Also, while viewing the Factory Modes, earlier encountered errors stored in memory can be printed.

To access the TESTS MENU, press the 5 Key from the MAIN menu. The TESTS MENU appears showing eight choices as shown in Figure 39.

NONACTIVE KEYS

The 1 and 4 key selections are nonactive and have no function assigned. Pressing either key will only result in a audible beep. (Refer to Figure 39).

Figure 39



Entering the Tests Menu

DICKEY-john®

VIEWING COM1 PORT (2 KEY)

The COM1 TEST screen displays the configuration of the COM1 (communications) port. Values are set from the Setup Mode to match the driving device (i.e. printer) of the port.

To access the COM1 port:

- Access the COM1 Menu by pressing the 2 key from the TESTS MENU (Refer to Figure 40). The COM1 TEST screen will appear displaying the COM1 Serial Data Setup status (BAUD RATE, DATA BITS, PARITY, STOP BITS) on the left half of the screen. The COM1 Control Lead Setup status displays on the right side displaying the Buffer Status. The BUFFER line indicates either BUSY or EMPTY.
- 2. Verify that a printer is connected and press the Print key to print a sample of the entire TRANSMIT CHARACTER SET. When completed, press the **Previous Page** key to return to the TESTS MENU.

Figure 40

Entering the Com1 Test Menu



VIEWING COM2 TRANSMIT PORT (3 KEY)

The COM2 TRANSMIT TEST screen displays configuration of the COM2 TRANSMIT port. Values are set from the Setup Mode to match the communicating device (i.e. computer) of the port.

To configure the COM2 TRANSMIT port.

1. From the TESTS MENU, select COM2 TRANSMIT TEST using the3 key (Refer to Figure 41). The screen displays the COM2 Mode, Serial Data Setup (BAUD RATE, DATA BITS, PARITY, STOP BITS), Buffer Status and the COM2 Control Lead Setup. If the buffer contains information, the line will indicate BUSY, otherwise EMPTY is indicated. Pressing the **Print** key produces a printout of the character set through the COM2 port. Press the **Previous Page** key to return to the TESTS MENU.



Print key



Previous Page key



Entering the Com2 Transmit Test Menu



VIEWING DISPLAY CHARACTERS (5 KEY)



Previous Page key

Typical characters are used to identify operator generated messages (Refer to Figure 42).

To access displayable characters:

1. From the TEST MENU, press the 5 key to select DISPLAY that shows available characters, symbols, and digits. After viewing, press the **Previous Page** key to return to the TESTS MENU.

Figure 42

Display Screen



DICKEY-john*



The Diagnostic Mode accesses unit operating values. Most of these readout values relate to the test cell circuitry.

To access Diagnostic details:

- 1. Press the 6 key to enter the DIAGNOSTIC MODE (Page 1). Two pages of information are available as shown in Figure 43.
- 2. Press the number key next to each selection, one key at a time. Each value calculates and appears next to the line title. If the printer is connected, a printout of each result is automatically printed (Refer to Figure 44 for page 1 printout and Figure 45 for page 2 printout).

The meaning of each line is as follows:

- 1: Vnull null voltage from cell board.
- 2: Vtemp temperature diode voltage from cell board.
- 3: Cwgt counts from weight oscillator.
- 4: Ce/f ratio of empty cell weight oscillator counts to full cell weight oscillator counts. When 4 key is pressed, item 3 blanks to recalculate. A message states PLEASE WAIT followed by PLACE WEIGHT ON CELL THEN PRESS ENTER. A test weight is placed on top of the test cell and then the Enter key is pressed. After a short wait, the results display.
- 5: SCD1 D1 (conductivity) self-check 379 ±20. When either the 5 or 6 key is pressed, these two items calculate.
- 6: SCD2 D2 (capacitance) self-check- 2500 ±20.
- 7: Temp Computed cell temperature in degrees C.

Figure 43

Diagnostic Mode Screens





Next Page key

3. Press the **Next Page** key to view Page 2 of the DIAGNOSTIC MODE. The display screen is divided into two columns; 5 kHz and 2 MHz (See Figure 43).



Enter key





Previous Page key

5 kHz (D1) – Left half of screen

1: Vc - Voltage of 5 kHz oscillator.

2: Vc-Vn - Voltage of 5 kHz oscillator minus null voltage.

3: Vsc- Voltage of 5 kHz oscillator with on-board test

network engaged.

4: Vsc - Vn - Voltage of 5 kHz oscillator with the on-board test network engaged minus null voltage.

2 MHz (D2) – Right half of screen

5: Vc - Voltage of 2 MHz oscillator.

6: Vc - Vn - Voltage of 2 MHz oscillator minus null voltage.

7: Vsc - Voltage of 2 MHz oscillator with on-board test network engaged.

8: Vsc - Vn - Voltage of 2 MHz oscillator with on-board test network engaged minus null voltage.

When finished, press the **Previous Page** key two times to return to the TESTS MENU.

Figure 44

Diagnostic Mode Page 1 Printout

VNULL: VTEMP: CWGT: CE/F: S/N: SCD1: SCD2: NEC: KEC: KEC:	0,0496 1,1838 38908 1,05377 0261-2360 378 2499 0,04941 1,50218 1,22890
SCD2: NEC: KEC: KSC: MEC: MSC:	2499 0,04941 1,50218 1,22889 1,40597 0,43373
TEMP:	18,26
\sim	\sim

Figure 45







USING TEST NETWORK MODE (7 KEY)

This mode is used by trained authorized personnel in troubleshooting. The screens and printout are shown for reference only.

To Access the TEST NETWORK MODE:

- 1. Press the 7 key TEST NETWORK MODE (See Figure 46). The message PLEASE WAIT will appear followed by the PRESS LOAD TO BEGIN.
- 2. Press the **Load** key to start the measuring cycle. When the test completes, the results display and automatically print if a printer is connected (See Figure 47).



Load key







Figure 47

Test Network Mode Printout

DICKEY-JOHN, CORP. GAC2100 5/18/96 09:18:59 ID:		
5/18/96 09:18:59 ID:		DICKEY-JOHN, ORP. GAC2100
ID:	5/18/96	09:18:59
NE: 400	ID:	
NI 430 NE: 490 D1 0,0 D1CF: 15000 D2: 14,4 D2CF: 14853 D3: 0,0 D4: 204,3 D4TF: 20,43 D4TE: 20,39 SELECT TEST NUMBER	NF: NE: D1 D1CF: D2CF: D2CF: D3: D4: D4TF: D4TE: SELEC	490 490 0,0 15000 14,4 14853 14853 0,0 204,3 20,43 20,39 T TEST NUMBER

DICKEY-john®



Print key



Previous Page key

VIEWING/PRINTING PARAMETERS (8 KEY)

The PARAMETERS MENU (8 key) operates similar to the SERVICE MENU, located in the SETUP MENU. The major difference is the inability to edit values on these screens. However, the VIEW FACTORY MODE screen is useful in viewing the number of currently stored errors in memory.

To view the PARAMETERS MENU:

- 1. Select the PARAMETERS MENU by pressing the 8 key from the TESTS MENU (See Figure 48). This menu allows the operator to view the unit parameters.
- 2. To print the values, select any menu and press the **Print** key. A printout of all the screens results (See Figure 49). Item 6 also prints all stored errors.
- 3. When finished viewing any screen, press the **Previous Page** key to return to the TESTS MENU.



Viewing Parameters Menu





Parameter Values Printout

\frown	
5/18/96	09:59:46
S/N: 0261-2	23601
REF TEMP REF VOLT:	: 25,00 1,1560
EC COUNT	: 40399
WEIGHT TA 1 0,0000 2 61,500 3 120,400 4 181,800 5 240,600 6 300,300	ABLE: 1,00000 1,02754 1,05453 1,08340 1,11179 1,14108
CAC: 1,000	0
IC: 1122	
BURN IN M	ODE: OFF
ERROR ST	ORING: ON
CURRENT ERROR 2 05/04/96 11 05/17/96 10 5 50134	ERRORS 8:42:2 6:50:4
SELECT TE	EST NUMBER





Load key





Previous Page key

OPERATOR MAINTENANCE

CELL CLEANING - DAILY

- 1. From the MAIN menu press the 5-key to display the TESTS menu. Press the 6-key to display the DIAGNOSTICS menu.
- 2. At the DIAGNOSTICS menu, press the **Load** key. The hopper doors open and remain open until **Unload** is pressed.
- 3. Use the brush supplied in the accessory kit to clean the inside of the grain cell.
- 4. Slide the GAC[®] 2100 G/GI to the front of the worktable and remove the grain drawer. Rach up inside and very carefully clean the strike off arm spring.

ACAUTION

Before performing Step 5, ensure that all fingers or other foreign objects are not in the cell area.

 Press the Unload key. The cell will rotate 180 degrees to dump the residue from the cell cleaning procedure. The hopper doors close. Pressing the Previous Page key twice returns the display to the MAIN menu.

FUSE REPLACEMENT

- 1. Disconnect the line cord from the rear of the unit. Using a small blade screwdriver or similar tool, remove the cover and fuse block assembly.
- Replace the blown fuse with a 1 Amp Slo-Blo fuse (Fuse Size = 5mm (0.197") dia. x 20 mm (0.787") long.
- 3. Reinstall fuse block assembly and cover.

Figure 50





DICKEY-john®

ERROR CODES

The GAC[®] 2100 G/GI contains a microprocessor to control the grain sample measurements and self-checks to determine the integrity of the internal electronics. If any limit is exceeded or malfunctions occur, an error code displays.

ERROR 1

Probable Cause:

The empty cell self check reading is out of limits.

Corrective Action:

Dump the cell and ensure it is empty.

ERROR 2

Probable Cause:

The empty cell weight reading is out of limits.

Corrective Action:

Dump the cell and ensure it is empty.

ERROR 3

Probable Cause:

The unload mechanism has malfunctioned.

Corrective Action:

Check the cell position and ensure it is free to operate.

ERROR 4

Probable Cause:

The strike-off mechanism has malfunctioned.

Corrective Action:

Ensure the cell and strike-off arm are free to operate.

ERROR 5

Probable Cause:

An invalid calibration is referenced.

Corrective Action:

Reselect the same grain or select a new grain and try again. If the condition continues, contact DICKEY-john Service Department.



ERROR 6

Probable Cause:

The moisture reading exceeds the upper limit.

Corrective Action:

Select a new grain and try again.

ERROR 7

Probable Cause:

The moisture reading exceeds the lower limit.

Corrective Action:

Select a new grain and try again.

ERROR 8

Probable Cause:

The empty cell ambient temperature has exceeded its limit.

Corrective Action:

Allow the GAC[®] 2100 G/GI to reach normal operating temperature.

ERROR 9

Probable Cause:

The sample temperature exceeds its limit.

Corrective Action:

Allow the sample to either cool down or warm up, then retry.

ERROR 10

Probable Cause:

The difference in temperature between empty cell ambient and full cell condition exceeds its limits.

Corrective Action:

Change the appropriate conditions and try again.

ERROR 11

Probable Cause:

Test weight exceeds stored average bulk density value by more than + 20 percent or -30 percent.

Corrective Action:

Use a more representative sample.

ERROR 12

Probable Cause:

A system memory error occurs.

Corrective Action:

Turn the power off for at least five seconds and then back on again. If the condition continues, contact DICKEY-john Service Department.

ERROR 13

Probable Cause:

An on-board RAM memory error has occurred.

Corrective Action:

Turn the power off for at least five seconds and then back on again. If the condition continues contact the supervisor.

ERROR 14 - NOT USED

ERROR 15

Probable Cause:

The real time clock battery is faulty.

Corrective Action:

Contact DICKEY-john Service Department.





APPENDIX A COMMUNICATION PORTS

COM1 PORT OPTIONS

The COM1 Port is a serial port conforming to the EIA/TIA RS-232-E standard and configured as Data Communication Equipment (DCE). Data out is RS-232-E levels and parity with adjustable format, data setup and control.

Table 1. COM1 Programmable Options		
Function	Option	Hex
Output Format	PRM 20	
	PRN 80	
	TST 20	
	PRN 81	
	CSV 01	
	PRN 21	
Line Terminator	LF (Line Feed)	0A
	CR (Carriage Return)	0D
	LF & CR	0A, 0D
	CR & LF	0D, 0A
Data Setup and Baud Rate	300	
	600	
	1200	
	2400	
	4800	
	9600	
	19200	
Control Busy	Ignore	
	Active (+)	
	Active (-)	
Control DTR	Ignore	
	Active	

COM1 PIN ASSIGNMENT DETAILS

COM1 is configured as Data Communication Equipment (DCE).

- **Protective Ground (pin 1)** electrically connected to frame.
- Received Data (pin 3) transmitted data from the GAC[®] 2100 G/GI.
- Data Set Ready (pin 6) ON (positive voltage) whenever the GAC[®] 2100 G/GI is powered.
- Signal Ground (pin 7) common reference for all other COM1 signals.
- Received Line Signal Detector (pin 8) ON (positive voltage) whenever the GAC[®] 2100 G/GI is powered.



- Busy (pin 11) gates the transmission of data from the GAC[®] 2100 G/GI. The COM1 CONTROL LEAD SETUP can be configured ACTIVE (+), ACTIVE (–) or IGNORE. When the line is inactive, up to 2048 characters are buffered. If not connected, configure this line as IGNORE.
- Data Terminal Ready (pin 20) indicates the connected device is not available. The COM1 CONTROL LEAD SETUP screen configures this line as ACTIVE or IGNORE. If configured ACTIVE, an OFF condition (negative or zero voltage) inhibits the transmission of data from the GAC[®] 2100 G/GI. Results display but not buffered for output on COM1. If not connected, this line should be configured as IGNORE.

Table 2. COM1 Connector Pins Identification			
Pin Number	Designation	Signal Name	Direction
1	AA	Protective Ground	NA
3	BB (RxD)	Received Data	To Device
6	CC (DSR)	Data Set Ready	To Device
7	AB	Signal Ground	NA
		Received Line Signal	
8	CF (RLSD)	Detector	To Device
11		Busy	To GAC2100
20	CD (DTR)	Data Terminal Ready	To GAC2100



COM2 PORT OPTIONS

Table 3. COM2 Programmable Options		
Function	Option	Hex
Output Format	PRM 20	
	PRN 80	
	TST 20	
	PRM 81	
	CSV 01	
	PRN 21	
Line Terminator	LF (Line Feed)	0A
	CR (Carriage Return)	0D
	LF & CR	0A, 0D
	CR & LF	0D, 0A
Data Setup	300	
Baud Rate	600	
	1200	
	2400	
	4800	
	9600	
	19200	
Data Bits	7 or 8	
Parity	None, Even, or Odd	
Stop Bits	1 or 2	
Control MODE	Transmit Only	
Control RTS	On (+)	
	Auto	
Control DRS	Ignore	
	Active	

COM2 PIN ASSIGNMENT DETAILS

COM2 is configured as Data Terminal Equipment (DTE).

Protective Ground (pin 1) – electrically connected to frame. Transmitted Data (pin 2) – data is transmitted from this pin. Received Data (pin 3) – data is received on this pin.

Request To Send (pin 4) – configured ON (positive voltage) whenever unit is powered or to be ON (positive voltage) only when there is data in the COM2 buffer to be transmitted. Use the COM2 CONTROL LEAD SETUP screen to configure this line.

Clear To Send (pin 5) – gates transmission of data from the $GAC^{\ensuremath{\mathbb{R}}}$ 2100 G/GI. This line must be ON (positive voltage) to obtain data on COM2. If an appropriate signal is not available from the connected device, connect this pin directly to Request To Send (pin 4). When the line is inactive, up to 2048 characters are buffered.


Data Set Ready (pin 6) – notifies the GAC[®] 2100 G/GI that the connected device is not available. The COM2 CONTROL LEAD SETUP screen configures this line as active or ignored. If configured as active, an OFF condition (negative or zero voltage) inhibits the transmission of data from the GAC[®] 2100 G/GI. Results display but not buffered for output on COM2. If not connected, this line should be configured as ignored. **Signal Ground (pin 7)** – common ground reference for all other COM2 signals.

Received Line Signal Detector (pin 8) – currently not used. Data Terminal Ready (pin 20) – ON (positive voltage) whenever the GAC[®] 2100 G/GI is powered. Ring Indicator (pin 22) – currently not used.

Table 4. COM1 Connector Pins Identification			
Pin Number	Designation	Signal Name	Direction
1	AA	Protective Ground	NA
2	BA (TxD)	Transmitted Data	To Device
3	BB (RxD)	Received Data	To GAC2100
4	CA (RTS)	Request to Send	To Device
5	CB (CTS)	Clear to Send	To GAC2100
6	CC (DSR)	Data Set Ready	To GAC2100
7	AB	Signal Ground	NA
		Received Line Signal	
8	CF (RLSD)	Detector	To GAC2100
20	CD (DTR)	Data Terminal Ready	To Device
22	CE (R1)	Ring Indicator	To GAC2100

Dealers have the responsibility of calling to the attention of their customers the following warranty prior to acceptance of an order from their customer for any DICKEY-john product.

DICKEY-john[®] WARRANTY

DICKEY-john warrants to the original purchaser for use that, if any part of the product proves to be defective in material or workmanship within one year from date of original installation, and is returned to DICKEY-john within 30 days after such defect is discovered, DICKEY-john will (at our option) either replace or repair said part. This warranty does not apply to damage resulting from misuse, neglect, accident, or improper installation or maintenance. Said part will not be considered defective if it substantially fulfills the performance expectations. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES OF MERCHANTABILITY, FITNESS FOR PURPOSE, AND OF ANY OTHER TYPE, WHETHER EXPRESS OR IMPLIED. DICKEY-john neither assumes nor authorizes anyone to assume for it any other obligation or liability in connection with said part and will not be liable for consequential damages. Purchaser accepts these terms and warranty limitations unless the product is returned within fifteen days for full refund of purchase price.

For DICKEY-john Service Department, call 1-800-637-3302 in either the U.S.A. or Canada



Headquarters:

5200 Dickey-john Road, Auburn, IL 62615 TEL: 217-438-3371, FAX: 217-438-6012, WEB: www.dickey-john.com

Europe:

DICKEY-john Europe S.A., 165, boulevard de Valmy, 92706 - Columbes - France TEL: 33 (0) 1 41 19 21 80, FAX: 33 (0) 1 47 86 00 07